

GET WRAPPED UP IN AUTUMN SAVINGS!

Flexible Voucher Option! Secure your budget now and use your voucher to run a project within the next 12 months.

1-Step 16S and ITS Microbiome analysis promotion and voucher option

1-Step Bacterial microbiome analysis or 1-Step Fungal microbiome analysis:

16S V4 region, Primer combination: 515F – 806R¹

16S V4-V5 region, Primer combination: 515F – 909R¹

16S V3-V4 region, Primer combination: 341F – 806R⁴

ITS ITS 1 region, Primer combination: ITS1F – ITS2²

Price example:

Package A: 192 samples only 31,20 € per sample and one region

Package B: 96 samples only 46,25 € per sample and one region

Package C: 48 samples only 74,50 € per sample and one region

Package D: 24 samples only 135,00 € per sample and one region

Package E: 10 samples only 176,50 € per sample* and one region

(*limited to V4 region)

Eukaryotes (microbial) analysis:

V8-V9 region, Primer combination: 18S-1422f (5'-ataacaggtctgtgatgcctt-3')
– 18S_1797R (5'-ccttcygcaggttcacctac-3')⁵

V9 region, Primer combination: Illumina_Euk_1391f (5'-GTACACACCGCCCGTC-3')
– Illumina_EukBr_1510r (5'-TGATCCTTCTGCAGGTTACCTAC-3'), blocking primer optional³;

V9 region samples could not be combined with ITS, V8-V9 or 16S samples and need to be sequenced in a separate MiSeq run because of smaller amplicon size!

Send us your gDNA samples and make use of our "All in One" service:

- Quality control
- Single step amplicon generation with reduced bias
- Double indexing, quality check, quantification, normalization and pooling of amplicons
- Illumina MiSeq / NextSeq 2000™ sequencing package A-D: 2 x 300 nt paired-end sequencing with V3 chemistry
- Illumina MiSeq sequencing package E: 2 x 250 nt paired-end sequencing with V2 chemistry nano flow cell
- Output package A-C: 20-30 M reads (including 25 % PhiX to balance the composition of bases)
- Output package D: 750 K reads (including 25 % PhiX to balance the composition of bases)
- De-multiplexing of reads
- Data delivery via FTP server

^{1,2,3,4,5} Primer sequence references: Please see page 4!

StarSEQ

Prices without tax. Cannot be combined with other promotions or discounts. Valid until 31.12.2023.

GET WRAPPED UP IN AUTUMN SAVINGS!

Flexible Voucher Option! Secure your budget now and use your voucher to run a project within the next 12 months.

2-Step 16S and ITS Microbiome analysis

2-Step Bacterial and Archaea Microbiome analysis:

16S, V1-V2 region, Primer combination 27F (5'-AGAGTTTGATCMTGGCTCAG-3') – 338R (5'-TGCTGCCTCCCGTAGGAGT-3')

2-Step Fungal (18S) microbiome analysis:

ITS, ITS 1/2 region, Primer combination: ITS-u2 (5'-GAAYCATCGARTCTTTGAACGC-3') – ITS-p4 (5'-CCGCTTAKTGATATGCTTAAA-3')

ITS, ITS 2 region, Primer combination: ITS3F (5'-GCATCGATGAAGAACGCAGC-3') – ITS4R (5'-TCCTCCGCTTATTGATATGC-3')

2-Step Barcoding and Custom Amplicon analysis

2-Step Barcoding:

cytochrome b, cytb1 (5'-ccatccaacatctcagcatgatgaaa-3') – cytb2 (5'-gccctcagaatgatattgtcctca-3')

cytochrome oxidase I, CO I E (5'-ccagagattagagggaatcagt-3') – CO 1 F (5'-cctgcaggaggaggagaycc-3')

or 2-Step Custom Amplicons:

design of custom primers with Illumina specific overhangs

Request for offer: microbiome@starseq.com

Optional: gold standard bioinformatic analysis of 16S/18S/ITS sequences (QIIME2)

The bioinformatics analysis pipeline consists of quality control, pre-processing of reads, taxonomy identification and visualization, calculation of alpha and beta-diversity metrics.

Request for offer: microbiome@starseq.com



GET WRAPPED UP IN AUTUMN SAVINGS!

Flexible Voucher Option! Secure your budget now and use your voucher to run a project within the next 12 months.

16S and ITS Microbiome analysis



Sequencing of your ready to load 16S/ITS/18S libraries

- Illumina MiSeq / NextSeq 2000™ sequencing: 2 x 300 nt paired-end sequencing with V3 chemistry
- Output 20-30 M reads (including 25 % PhiX to balance the composition of bases)
- De-multiplexing of reads
- Free 16S metagenomics analysis (Illumina App)
- Data delivery via FTP server

Only 3250 €

Optimize your research with our 16S/ITS workflow solutions



From sample prep to library prep to sequencing and bioinformatics. StarSEQ offers a complete and comprehensive workflow for various 16S/ITS applications. Sequencing of the hypervariable regions of the 16S rRNA gene is the gold-standard for identification and classification of bacterial communities within a given sample.

StarSEQ is using established gold-standard primer sequences to target regions covering the hypervariable regions V3-V4, V4 and V4-V5. In contrast to common used methods, StarSEQ has designed its own single-step technique for simultaneous amplification, indexing and adapting of samples to reduce the bias caused during PCR.

StarSEQ

Prices without tax. Cannot be combined with other promotions or discounts. Valid until 31.12.2023.

GET WRAPPED UP IN AUTUMN SAVINGS!

**Flexible Voucher Option! Secure your budget now
and use your voucher to run a project within the next 12 months.**

16S and ITS Microbiome analysis

Applications:

- Environmental Metagenomics (soil, water, air, biofilms, complex organic communities)
- Human or Animal Microbiome (skin, stool, gut, blood, swab)
- Monitoring of animal health (e.g. alternatively or supplementarily to FELASA-Test)
- Sterility Monitoring
- Detection of Contamination
- Monitoring of Biogas Plant
- Biosafety Monitoring
- Food Quality
- Clinical Samples

**You have no time/staff for routine preparation
of gDNA from your collected samples?**

**As an additional service we offer also the
purification of gDNA starting from 18 €.**

Primer sequence references:

¹ Apprill A, McNally S, Parsons R, Weber L. 2015. Minor revision to V4 region SSU rRNA 806R gene primer greatly increases detection of SAR11 bacterioplankton. *Aquat Microb Ecol* 75:129–137. Parada AE, Needham DM, Fuhrman JA. Every base matters: assessing small subunit rRNA primers for marine microbiomes with mock communities, time series and global field samples. *Environ Microbiol*. 2016;18: 1403-1414.

¹ Parada, A. E., Needham, D. M., & Fuhrman, J. A. (2016). Every base matters: assessing small subunit rRNA primers for marine microbiomes with mock communities, time series and global field samples. *Environmental Microbiology*, 18(5), 1403–1414. <https://doi.org/10.1111/1462-2920.13023>

¹ Walters, W., Hyde, E. R., Berg-Lyons, D., Ackermann, G., Humphrey, G., Parada, A., ... Knight, R. (2016). Improved Bacterial 16S rRNA Gene (V4 and V4-5) and Fungal Internal Transcribed Spacer Marker Gene Primers for Microbial Community Surveys. *mSystems*, 1(1), e00009-15. <https://doi.org/10.1128/mSystems.00009-15>

² White, T. J., T. Bruns, S. Lee, and J. W. Taylor. 1990. Amplification and direct sequencing of fungal ribosomal RNA genes for phylogenetics. Pp. 315-322 In: PCR Protocols: A Guide to Methods and Applications. Academic Press, New York, NY. Gardes, M., and T. D. Bruns. 1993. ITS primers with enhanced specificity for basidiomycetes - application to the identification of mycorrhizae and rusts. *Mol. Ecol.* 2: 113-118.

³ Amaral-Zettler, L. A., McCliment, E. A., Ducklow, H. W., & Huse, S. M. (2009). A Method for Studying Protistan Diversity Using Massively Parallel Sequencing of V9 Hypervariable Regions of Small-Subunit Ribosomal RNA Genes. *PLOS ONE*, 4(7), e6372. Retrieved from <https://doi.org/10.1371/journal.pone.0006372>

⁴ Klindworth A., Pruesse E., Schweer T., Peplies J., Quast C., Horn M., et al. . (2013). Evaluation of general 16S ribosomal RNA gene PCR primers for classical and next-generation sequencing-based diversity studies. *Nucleic Acids Res.* 41, 1–11. 10.1093/nar/gks808
see also https://support.illumina.com/downloads/16s_metagenomic_sequencing_library_preparation.html

⁵ Bradley, I. M., Pinto, A. J., & Guest, J. S. (2016). Design and Evaluation of Illumina MiSeq-Compatible, 18S rRNA Gene-Specific Primers for Improved Characterization of Mixed Phototrophic Communities. *Applied and Environmental Microbiology*, 82(19), 5878 LP-5891. <https://doi.org/10.1128/AEM.01630-16>
see also https://support.illumina.com/downloads/16s_metagenomic_sequencing_library_preparation.html



Prices without tax. Cannot be combined with other promotions or discounts. Valid until 31.12.2023.